

# Defoam 2000

## Defoaming agent on fluorosilicone basis for the distilling of perc and hydrocarbon solvents

---

- Properties:**
- destroys foaming reliably and prevents „boiling over“ of the distilling vessel
  - acts successfully also in silicone containing liquors
  - prevents foaming effectively
  - does not contain any volatile components or any agents harmful to the solvent system

**Application:** Give **Defoam 2000** into distilling vessel before beginning of distillation. Since the active substance of **Defoam 2000** remains in the distilling vessel the defoaming effect continues as long as the distilling vessel is cleaned.

Dosing recommendation

*Low foam formation*

10 - 20 ml per 100 l solvent (12.5 - 25 fl oz/10 gal)

*Strong foam formation*

20 - 40 ml per 100 l solvent (25 - 50 fl oz/10 gal)

**Technical Data:**

Density 20 °C (68 °F)	0.87 g/ml
-----------------------	-----------

**Hints:**

Storage

Store container tightly closed and upright. The product can be stored for at least 24 months in its original sealed packing.

1104